

WVU IACUC Model Guidance Sheet: Use of Streptozotocin (STZ) for Diabetes Induction in Animals

Purpose

STZ is a commonly used chemical agent for the induction of diabetes in research animals. This compound poses human health and safety risks and precautions ***must*** be taken. In addition, diabetes model is a chronic disease model which monitoring, supportive care and humane endpoints ***must*** be considered to ensure that appropriate animal welfare is maintained.

Definitions

Streptozotocin (STZ)- A compound derived from *Streptomyces achromogenes* with antineoplastic properties which target the insulin producing beta cells of the pancreas in mammals.

Carcinogen- A substance capable of causing cancer.

Teratogen- A substance which can cause malformation of an embryo/fetus.

Cytotoxic- A substance which is toxic to cells.

Guidance

1. Introduction
 - a. STZ is a highly toxic compound. It is considered cytotoxic to pancreatic beta cells, suspected carcinogen, mutagen, teratogen, and reproductive toxin. There is no established safe occupational exposure limit.
 - b. Potential routes of exposure include inhalation, ingestion, accidental injection, and dermal absorption.
 - c. All safety procedures outlined by Environmental Health and Safety ***must*** be adhered to when working directly with this compound or handling animals which received this compound. Please follow the IACUC SOP "*Hazardous Chemicals Used with Animals*"
 - d. Animals are considered to be hazardous for 72 hours after last administration of STZ. After the 72-hour period, animals can be transferred back into standard housing rooms.
2. Reagents
 - a. STZ is commonly obtained as a non-pharmaceutical grade compound (*see IACUC Policy "Non-Pharmaceutical-Grade Substances Used in Animals"*) from vendors such as Sigma-Aldrich®.
 - b. STZ is typically mixed in a citrate buffer solution to obtain a final compound with a pH of 4.5 - 5.5.
 - c. Compounds ***must*** be prepared using sterile technique.
 - d. All reagents used and a description of the compounding process ***must*** be described in the animal use protocol.
 - e. Compounding considerations:
 - i. STZ should be stored at -20° C.

- ii. STZ container should be protected from light due to potential degradation.
 - iii. STZ is unstable and should be mixed with citrate buffer immediately prior to injection. It ***must*** be prepared fresh each time and given within 5 minutes of mixing.
3. Common Dosing Regimens
- a. STZ can be used to induce either Type I or Type II diabetes depending on the administration regimen used.
 - b. Sensitivity to compound can vary with animal age, strain/stock, and sex used. A pilot study may be necessary to determine appropriate treatment regimen. Below are some commonly used protocols to provide a starting point, but regimens may need to be adjusted based on response.
 - c. Some diabetes induction protocols require a short fasting period prior to STZ administration.
 - i. Standard fasting is 4 hours for mice and 6-8 hours for rats.
 - d. Rats
 - i. Type I diabetes model
 - Single dose 40-70 mg/kg (65 mg/kg commonly used dose) in rats 8-10 weeks of age
 - ii. Type II diabetes model (STZ + Nicotinamide)
 - Nicotinamide (230 mg/kg IP) followed by STZ (65 mg/kg IV)
 - Nicotinamide protects against complete beta cell toxicity
 - iii. Type II diabetes (fat fed STZ model)
 - High fat diet (60% fat by caloric content) for 3 weeks prior to STZ administration
 - STZ given on day 22 (40 mg/kg IP)
 - e. Mice
 - i. Repeated low dose Type I diabetes model
 - 40 mg/kg IP once daily for 5 consecutive days
 - This model has fewer toxic effects
 - ii. Single high dose Type I diabetes model
 - Single IP injection 200 mg/kg
4. Prior to administration of STZ
- a. Ensure animals are transferred into a chemical hazard room and Office of Laboratory Animal Resources (OLAR) staff has appropriate documentation (Animal Hazard Form and MSDS). Animals receiving hazardous chemicals ***must*** be housed in a designated chemical hazard room and cannot remain in the standard housing. Contact OLAR husbandry managers for guidance on this process.
 - b. STZ injections and all animal manipulations (during first 72 hours post-STZ administration) ***must*** be performed under a certified class 2A or 2B biosafety cabinet or chemical fume hood.
 - c. Animals ***must*** be properly restrained for injection. Chemical restraint reduces risk of potential exposure of personnel to hazardous agents.
 - d. All animals should be weighed on day of experiment to ensure accurate dosing of STZ compound.
5. STZ administration
- a. STZ should be administered as described in the approved IACUC animal use protocol.
 - b. After injection, animals should be returned to their cage with adequate food and 10% sucrose water (prevents post-injection hypoglycemia). Sucrose water should be maintained on cage for 48-72 hours post-injection.

6. Post-STZ administration monitoring and considerations
 - a. Animal use protocol should outline frequency and method used for blood glucose (BG) testing.
 - i. Fasting is typically performed prior to BG testing. The NIH and Animal Models of Diabetic Complications Consortium (AMDCC) recommend a morning fast starting at 7 am with an afternoon blood collection. Fasting recommendation is 6 hours for mice and 6-8 hours for rats.
 - b. Animal use protocol should outline expected BG concentration once animal is diabetic and any compounds administered to reduce BG.
 - c. Monitor for clinical signs of hypoglycemia during first 48 hours.
 - i. Lethargic
 - ii. Unresponsive
 - iii. Contact veterinary staff to assist with treatment
 - d. Recommend daily monitoring for the first week post injection, and weekly monitoring once stable.
 - e. Diabetic animals have an increased water intake need, please ensure adequate water supply is available.
 - f. Diabetic animals may require an increased cage change frequency due to increased urination.
 - g. Clinical signs of diabetic animals:
 - i. Polydipsia (increase water intake)
 - ii. Polyuria (increase urination)
 - iii. Abdominal distention
 - iv. Dehydration
 - v. Weight loss and body condition loss
 - vi. Unkempt appearance
 - h. Diabetic animals can benefit from supportive care including nutritional support (soft chow/diet gel), food on ground/dish on cage floor, fluid support, heat support, group housing, nesting material/huts. A special care form should be submitted to OLAR staff describing additional care which is needed.
 - i. Routine monitoring (weighing, body condition score (BCS), hydration checks, BG measurements) and humane early endpoints ***must*** be described in the animal use protocol.

References

<https://currentprotocols.onlinelibrary.wiley.com/doi/epdf/10.1002/cpz1.78>